

MammoReport^{Plus} - MAMMOMAT 3000 Nova

SP

Installation Instructions

Settings for system including the Fuji CR Console

© Siemens AG 2003

The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights, including rights created by patent grant or registration of a utility model or design, are reserved.

Chapter	Page	Revision
All	All	01

Document revision level

The document corresponds to the version/revision level effective at the time of system delivery. Revisions to hardcopy documentation are not automatically distributed.

Please contact your local Siemens office to order current revision levels.

Disclaimer

The installation and service of equipment described herein is to be performed by qualified personnel who are employed by Siemens or one of its affiliates or who are otherwise authorized by Siemens or one of its affiliates to provide such services.

Assemblers and other persons who are not employed by or otherwise directly affiliated with or authorized by Siemens or one of its affiliates are directed to contact one of the local offices of Siemens or one of its affiliates before attempting installation or service procedures.

	Page
1 Prerequisites	1 - 1
Purpose1 - 1
Important note.1 - 1
Validity.1 - 1
Components included1 - 1
2 Worksteps	2 - 1
Settings2 - 1
Final procedures.2 - 1
3 Installation protocol	3 - 1
Settings3 - 1
System environment.3 - 1
MAMMOMAT upgrade.3 - 1

This page intentionally left blank.

Purpose

This document describes extra settings for a system comprising the MAMMOMAT 3000 Nova, MammoReport^{Plus}, Fuji CR Console and FCR5000 MA reader.

Important note

After the instructions in this document have been performed, the installation protocol at the end of this document must be signed in by the CSE and a Fuji Technician to confirm the installation. A copy of the signed installation protocol shall be sent to:

*Siemens AG
Medical Solutions
SP WH - Product Quality Group
Henkestrasse 127
D - 91052 Erlangen*

Validity

The instructions in this document must only be performed for a system comprising the MAMMOMAT 3000 Nova, MammoReport^{Plus}, Fuji CR Console and FCR5000 MA reader. The worksteps have to be performed by a Fuji Technician.

Components included

The package contains the following parts:

Part no.	Name
66 53 997	Workflow instructions, Supplement to the instructions for use
66 54 003	This instruction

This page intentionally left blank.

Settings

The following settings shall be made for the Fuji CR Console in order to function properly with the other units in the system. The settings shall be performed by a Fuji Technician.

- The following six menus shall be setup:
 - CC left
 - CC right
 - MLO left
 - MLO right
 - ML left
 - ML right
- The six menus shall be based on the menus 0307 and 0317, but with the following image processing parameters:
 - PEM off
 - DRC on
 - TAS off
 - GA=1.8, GT=Q, GC=0.6, GS=0
 - MRB=G, MRT=P, MRE=0.4
 - EDR=AUTO

Final procedures

After the settings have been performed, fill in the Installation protocol at the end of this document and send it to:

*Siemens AG
Medical Solutions
SP WH - Product Quality Group
Henkestrasse 127
D - 91052 Erlangen*

This page intentionally left blank.

Settings

	Passed	Failed
The six menus have been setup:	<input type="checkbox"/>	<input type="checkbox"/>
The image processing parameters have been setup:	<input type="checkbox"/>	<input type="checkbox"/>

Date of installation:

Name and signature,
Fuji Technician:

System environment

	Passed	Failed
A Hardcopy printer or an Archive system (e.g. PACS) is connected to the system:	<input type="checkbox"/>	<input type="checkbox"/>

Name of Hardcopy printer or Archive system:

Name and signature, CSE:

MAMMOMAT upgrade

	Passed	Failed
The AEC tables of the MAMMOMAT are calibrated with the Upgrade kit AEC (Fuji) for MAMMOMAT:	<input type="checkbox"/>	<input type="checkbox"/>

Name and signature, CSE:

When this protocol has been completed, file a copy with the MammoReport^{Plus} system binder and send this protocol to:

*Siemens AG
Medical Solutions
SP WH - Product Quality Group
Henkestrasse 127
D - 91052 Erlangen*

This page intentionally left blank.